NOVACer[®] ZR

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*Instruction of Use
Zirconium Dioxide layering ceramic*

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Safety information: *Please read these instructions thoroughly before using the product. Improper use of the product and failure to observe the information provided may affect quality and reduce the lifetime of the dental prosthesis. Only trained staff (dentists or dental technicians) should be allowed to handle this product.*

– WARNING –

Applies to porcelains in powder form: Please open the packaging and use the product only with protective equipment. Wear safety goggles and a dust mask. Dust particles irritate the eyes, mucous membranes and skin, and injure the lungs.

Intended use: NOVACer® for zirconium dioxide ceramic (ZR) covers modular ceramic materials for veneering the zirconium framework.

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Frame:

The zirconium oxide copings are produced in the CAD/CAM procedure. Zirconium oxide not only has great strength, it also meets the highest standards of aesthetics and biocompatibility.

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*The labial and occlusal shape is build-up with
NOVACer® ZR dentine in the required shade.*

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*The dentin will be reduced and chromatic mamelon
material build-up.*

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*Control of the correct volume of the
dentin structure after firing.*

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*The missing areas will be build-up with incisal and
transparent porcelains.*

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Glaze firing after finishing.

The degree of glaze can be adjusted by polishing.

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* The firing parameters given above are guidelines, which must always be adjusted to suit the furnace used and the situation of the furnace. What is essential is getting the right firing result.

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*All tested materials conform to EN ISO 9693-1. The technical and
physical values quoted relate to samples produced in-house at our
own measuring instruments.*

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